



# HOT/COLD THERAPY

## More Healthy Life



**INTCO Protects Your Health**

**Jiangsu Intco Medical Products Co., Ltd.**

No.298 Yandunshan Road, Dagang Zhenjiang, Jiangsu, 212132 China

[www.intcomedical.com](http://www.intcomedical.com) [info@intco.com](mailto:info@intco.com) +86 400-050-6868 +86 0511-83174088



V 2024.1.0



**INTCO PROTECTS  
YOUR HEALTH**



01.INTCO INTRODUCTION

Company Profile ..... 03

Development History ..... 05

Certifications of INTCO ..... 07

02.INSTANT/DISPOSABLE COLD OR HOT PACK

Cold Pack ..... 09

Hot Pack ..... 10

03.REUSABLE HOT & COLD PRODUCTS

Gel Pack ..... 12

Gel Beads ..... 17

Clay Pack ..... 18

Massage Ball ..... 19

04.ICE BOX & ICE PACK

Ice Box ..... 21

Ice Pack ..... 22

05.SPA SERIES

Eye Mask & Face Mask ..... 24

Moist Heat Wrap ..... 25

Magic Warmer ..... 26

06.DISPOSABLE HEAT PATCH

Heat Patch ..... 28

07.HOME CARE SERIES

Cooling Patch ..... 30

Ice Bag ..... 31

Hot Water Bag ..... 32

Cool Gel Mat ..... 33

Cooling Gel Pillow ..... 34

08.MOM & BABY SERIES

Breast Gel Pack ..... 36

Perineal Instant Cold Pack ..... 37

Infant Heel Warmers ..... 38

09.PROMOTION & GIFT SERIES

Cartoon Hot/Cold Toy Pack ..... 40





# 01 INTCO INTRODUCTION

## Jiangsu INTCO

Jiangsu INTCO Medical Products Co., Ltd. was established in 2004 and located in Dagang District, Zhenjiang City, Jiangsu Province. It covers an area of 64.2 acres and has a total investment of 121 million US dollars. It is a subsidiary of INTCO Medical Technology Co., Ltd. (Stock Code: 300677).

The company's products have covered many fields such as hot and cold physical therapy care, sanitation and disinfection, medical health, oral care, food and drug cold chain packaging, etc., we also have FDA/CE/ISO13485/BSCI and other quality system and product system certificates, besides 12 kinds of new practical product patents.



PRODUCING LINE



FRONTDESK




SHOWROOM



WAREHOUSE



# DEVELOPMENT HISTORY



## 1993

- Established INTCO Medical Industries Inc.

## 1998

- Established Shanghai Marketing Center.

## 2003

- Established ECG and Hot&Cold therapy manufacturing base.

## 2004

- Established wheelchair manufacturing base in Zhenjiang.

## 2009

- Established disposable gloves manufacturing base in Zibo.

## 2010

- Established disposable gloves manufacturing base in Qingzhou.

## 2013

- Established Canada Marketing Center.

## 2016

- Established Germany Marketing Center.

## 2017

- Listed on Shenzhen Stock Exchange.
- Established disposable gloves manufacturing base in Huaibei.

## 2019

- Established disposable gloves manufacturing base in Vietnam.

## 2020

- Established disposable gloves manufacturing base in Jiangxi.
- Established disposable gloves manufacturing base in Anqing.
- Established CHP project base in Anqing.

## 2021

- Established Japanese Marketing Center.

## 2022

- Established INTCO Medical Intelligent Science & Technology Park.
- Established Malaysia Marketing Center.





# CERTIFICATIONS OF INTCO




# 43 FACILITY PATENTS



# 02 INSTANT/DISPOSABLE COLD OR HOT PACK


Instant cold pack & Instant hot pack can provide fast and effective cold or hot temperature & therapy at once . They are widely used in following fields: hospital, home care, first aid, sports or other outdoor activities.






# COLD PACK

First Aid/Sports Support/Home Care.External application only.  
For relief of pain and swelling caused by sprains, strains, contusions, minor burns.  
This product can be customized in accordance with any sample size.



# HOT PACK

Provides warmth at low temperatures. Improve blood circulation.  
Relaxes the muscles and relieve arthritis.  
This product can be customized in accordance with any sample size.




## Characteristic

Single cold application use without the drippy mess of ice cubes.  
No refrigeration needs, convenient and portable.  
Primary usage for outdoor activity injuries and first aid.

## Usage

Shake contents to bottom of bag.  
Locate inner pouch. Squeeze firmly to rupture inner fluid pouch.  
Shake briefly to mix contents.  
Wrap it in a towel or a piece of cloth for cold compress, apply to affected area.

Item Number	Size	Weight	Material
C127153NB	12.7*15.3CM	105G	PE&NYLON
C153229NE	15.3*22.9CM	260G	PE&NYLON



## Characteristic

No heating required.  
Ideal for emergencies, work instantly.  
Convenient and portable.  
Magnesium sulfate formula is available.

## Usage

Shake contents to bottom of bag.  
Locate inner pouch. Squeeze firmly to rupture inner fluid pouch.  
Shake briefly to mix contents.  
Wrap it in a towel or a piece of cloth for cold compress, apply to affected area.

Item Number	Size	Weight	Material
H127153B	12.7*15.3CM	104G	PE&NYLON
H153229D	15.3*22.9CM	230G	PE&NYLON



# 03 REUSABLE HOT & COLD PRODUCTS

Designed for cold or hot application, this product is specifically contoured for body parts to serve as a flexible ice pack or for soothing heat. Cold can help reduce pain and swelling, while heat can alleviate muscle aches and stiffness.

## GEL PACK

### Instruction

For hot use, put the pack flat in microwave to heat, the heating time depends on the different sizes.  
For cold use, freeze the pack at least 2 hours in freezer before use.

### Cold therapy application

Reduce the swelling and relief sprains, strains, contusions, minor burns, toothaches, simple headaches.

### Heat therapy application

Relax tight muscles and restore flexibility in a wide variety of ailments such as:  
Muscle pain & soreness; Muscle spasms & cramps;  
Arthritis / Menstrual cramps; Chronic neck or back stiffness.

			
Item Number	HC-GEL 0507	HC-GEL 0710	HC-GEL 0510
Size	12.7*17.8CM	17*25CM	12.7*25.4CM
Weight	144G	375G	320G
Material	PE&NYLON	PVC	Fabric

			
Item Number	HC-GEL 0511	HC-GEL 0410	HC-GEL 1015
Size	12*29CM	12.5*28CM	25.4*38.1CM
Weight	350G	175G	800G
Material	PE&NYLON	PVC compound knitted fabric	PVC compound spunlace





GEL PACK


Characteristic

Remain soft and flexible after freezing.  
Different sizes and colors are available.  
Reusable with elastic wrap.

Usage

For heat therapy, microwave or boil the pack by instructions. For cold therapy, freeze the pack for at least one hour. Wear the wrap first, then put the pack on the wrap for use.

Item Number	WRA-R 2540	WRA-R 4060	WRA-R 5060	WRA-R 6060	WRA-R 7080
Size	R2540	R4060	R5060	R6060	R7080
Weight	92G	565G	795G	910G	1360G
Material	TPU	TPU	TPU	TPU	TPU



			
Item Number	WRA 0816	WRA 06010	WRA 0712
Size	40*20CM	16.5*25.9CM	15.3*35.6CM
Weight	850G	320G	400G
Material	PVC	PVC	Fabric
Wrap	Fleece/Neoprene fabrication	Fleece/Neoprene fabrication	Fabric

				
	Item Number	WRA 0721	WRA 1014	WRA 1319
	Size	17*53.4CM	25.4*34.5CM	32*47CM
	Weight	910G	1360G	2270G
	Material	TPU	TPU	TPU



NEW!

# NEW GEL COLD HOT PACK

The New Cold Hot Pack is designed for different part of body, which helps relieve pain, reduce inflammation,improves blood circulation,as a natural painkiller, cold hot pack can help you recover from muscle and joint injuries.

FOOT COVER

Pain Relief From

- Muscle and ankle joint injuries.
- Ankle Sprains
- Leg Spasms
- Leg Tightness

HEADWRAP

Pain Relief From

- TMJ
- Migraine
- Teeth Extractions

	<div>NEW</div>			<div>NEW</div>	
Item Number	WAR0410	HC-GEL0404	HC-GEL0409	WRP0531	HC-GEL0475
Size	12*27CM	10.1*12CM	10.2*24.5CM	13*80CM	10*19CM
Weight	40G	60G	140G	50G	170G
Material	Fleece	Fabric	Fabric	Fleece	Fabric



NEW!

# SOLID GEL PACK



Instruction

**Hot Therapy**—With just a few seconds in your microwave, your compression brace will be ready to be placed as a hot pad to any area that could use some soothing.

**Cold Therapy**—Place your ice pack sleeve in the freezer and enjoy its supreme comforting and anti-inflammatory action whenever pain strikes.

MIGRAINE RELIEF ICE HEAD WRAP

- Headache and Migraine Relief.
- Hot/ Cold Therapy.
- Snug Fit Provides Compression.
- Freezes Fast, Stays Cold.
- Dark Light-Blocking Eye Cover.

HOT & COLD COMPRESSION SLEEVE

- Effective Pain Relief.
- Comfortable & Slip-On Design.
- 360° Coverage Therapy.
- Hot / Cold Therapy.

	<div>NEW</div>	<div>NEW</div>
Item Number	HC-SLE630	HC-SLE480/HC-SLE450
Size	16*24CM	26*11CM/27*11.5CM
Weight	630G	480G/450G
Material	Sodium Polyacrylate, Glycerin, Water and Polyurethane Fibre	





# GEL BEADS

Gel beads pack have good performance for pain relief, provide long lasting cold and hot therapy, innovative gel beads and various shapes fit the body well and really eye-catching.



# CLAY PACK



### Instruction

**Hot Therapy**—Heat in microwave for specified time  
Heat for additional 5 seconds as needed until desired temperature is reached.

**Cold Therapy**—Freeze for at least 2 hours.


Item Number	HC-GB0326	HC-GB0407	HC-GB2038
Size	9.4*68CM	11.5*19CM	10*25CM
Weight	120G	160G	250G
Material	PVC	PVC	PVC
Item Number	HC-GB0511	HC-GB0411	HC-GB0617
Size	18.5*25CM	11.8*28.6CM	17.5*43CM
Weight	175G	160G	560G
Material	PVC	PVC	PVC
Item Number	HC-GB2894	HC-GB4124	HC-GB2819
Size	10.5*26.5CM	41*24CM	20.5*38.5CM
Weight	160G	500G	430G
Material	PVC+Flannelette	PVC+Flannelette	PVC+Flannelette



- It utilizes a synchronized blend of all-natural, ceramic clay and non-toxic oils to create the ideal thermal transfer source.
- It capitalizes on the efficiency of radiant energy to create deep penetrating, therapeutic pain relief.
- It is safe, easy to use, highly efficient in the treatment of a wide variety of pain discomforts such as arthritis, sports injuries, muscle pain , joint stiffness, menstrual cramps and headaches.
- The hot performance is more effective than gel.

Item Number	WRA-CLAY 0612	WRA-CLAY 0916
Size	15*30CM	24*41CM
Weight	670G	1400G
Material	PE&NYLON+WRAP	PE&NYLON+WRAP





# MASSAGE BALL

This speical massage ball filled with gel enables it to safely retain either heat or cold for up to certain minutes by simply placing in the freezer or hot water.  
Two different texture options can be chosen --- smooth or textured shape.



Massaging can help ease aches, pains and tension in your back,arms and legs.  
The aim of massaging is to relax and stretch the muscles, thus improving circulation and a free flow of oxygen, which ultimately benefits your sense of well-being.

### Instruction

Put into boiling water for heat use or put into freezer for at least two hours for cold use.

Item Number	MB-6	MBS-6	MBS-7
Size	Diameter 6CM	Diameter 6CM	Diameter 7CM
Weight	105G	165G	260G
Material	PVC	Steel/ABS	Steel/ABS



# 04 ICE BOX & ICE PACK

For transportation of fish, fruit, vegetable, blood and vaccine which need lower temperature to keep fresh.  
Put it into lunch bag(box) to keep food & drinks cool when you have outdoor picnics or camping.








ICE BOX

- Ideal for cool boxes, cool bags, and commercial catering and fishing.
- The ultimate reusable gel ice box.
- Tough. Rugged. Keeps icy-cool for longer.
- Non toxic commercial grade cooling gel.
- Size suitable for daily use.

Instruction

Put the ice box in the refrigerator for at least 3 hours before use, until completely frozen.

			
Item Number	BOX 350A	BOX 750A	BOX 1000A
Size	11*16.5*3CM	17*21*3CM	20*26*3CM
Weight	350G	750G	1000G
Material	HDPE	HDPE	HDPE

			
Item Number	BOX 530A	BOX 200A	BOX 160A
Size	9.8*5.2*21.3CM	12.2*17.8*1.2CM	13.3*12.7*1.4CM
Weight	530G	200G	160G
Material	HDPE	HDPE	HDPE



ICE PACK

- Long Lasting: The powerful cooling gel mixture-filled ice pack chill faster retains coldness much longer than ice which can keep your food and drinks cold for longer!
- Reusable & Safe: Easy to clean and reusable, which is great replacement for traditional ice cubes.
- Suitable for Many Occasions: Great for school lunch, camping, picnics, boating, fishing, or any need to keep cold cooler.



Instruction

Put the ice pack in the refrigerator for at least 3 hours before use, until completely frozen.

			
Item Number	PACK 126A	PACK 98A	PACK 500A
Size	13*13.5CM	12.7*12.7CM	13.5*21.5CM
Weight	126G	98G	500G
Material	PVC	PE&NYLON	PE&NYLON



# 05 SPA SERIES

Conforms to your body to deliver soothing relief from: migraines, sinus infections, cosmetic procedures and more Freeze or microwave. Reusable and non-toxic Latex, lead and BPA-free.






## EYE MASK & FACE MASK

Promote blood circulation in the face or eye, relieve fatigue.  
Relief pain and subsidence of a swelling.  
Get instant restoration and provide soft touch over the face, which is perfect for firs aid no matter at home or in the office.  
Also a great gift on any day.

### Instruction

Place eye mask into the refrigerator 1 hour before using. Wear it about 15 minutes.

			
Item Number	EYE-GB 0309	EYE-GB 0310	FACE-GB 0812
Size	22.6*7.6CM	24.5*7.5CM	19*30.8CM
Weight	96G	65G	230G
Material	PVC	PVC	PVC

			
Item Number	EYE-GEL 0208	EYE-GEL 0310	EYE-GEL 0314
Size	21*6CM	24.5*7.5CM	36*9.5CM
Weight	45G	75G	152G
Material	PVA	PVC	PVC





# MOIST HEAT WRAP


Moist Heat has been used for hundreds of years by many cultures in the treatment of arthritis, rheumatism, joint stiffness and muscle spasms. Moist Heat is preferred because of ability to penetrate deeply and more quickly at lower temperatures than dry heat.

### Instruction


For hot use, put the pack flat in microwave oven to heat, the heating time is depended on the different sizes.

Item Number	MH 250	MH 850A	MH 850B	MH 420
Size	21*10.3CM	18.5*39CM	17.6*55CM	50.8*10.2CM
Weight	250G	850G	850G	800G
Material	Crystal velvet/Fleece	Crystal velvet/Fleece	Crystal velvet/Fleece	Fabric

### STEP.1



### STEP.2



### STEP.3



# MAGIC WARMER



Just simply click and the instant heat can be provided. The magic warmer can be used for circulation boosting, deep penetrating heat therapy, sports therapy for muscles, joints strains and inflammation. There is a set of SPA series, which including shoulder, waist, foot, hand and square one.

### Instruction

Bend the metal disc in the liquid a few times to activate heating. Wrap soft cloth and apply pack to desired area.

Re-activation Method:  
Heat the water in a container until the temperature reaches 90-100°C. Immerse the activated warmer in hot water, for 10 to 15 minutes. Stop heating, remove the warmer from the water. Keep in a safe place away from children. When the warmer cool down, it can be used next time.

Item Number	MAG 0711	MAG 0816
Size	19.2*28.8CM	41*25CM
Weight	700G	840G
Material	PVC	PVC

Item Number	MAG-HEART	MAG-STAR	MAG-CIRCLE
Size	10.6*11.5CM	9.5*10CM	10*10CM
Weight	100G	53G	70G
Material	PVC	PVC	PVC



# 06 DISPOSABLE HEAT PATCH

Using heat patch is the most convenient and safest way to provide skin-touching continuous and stable warm. It is portable and has great effect in preventing from frostbite in winter. Just one single piece can resist full day's cold.



HEAT PATCH



### Main features

It's disposable and easy to provide relief of minor muscular and joint aches and pains associated with over exertion, strains, sprains and arthritis, period pain, etc.  
It will gradually get warm and last for at least 6/8 hours.

### Instruction

Open the outer pouch, leaving the internal pack to react with air. It needs some time to reach its therapeutic temperature. Please refer to the instruction for use for specific product usage.

Item Number	L095130A16	L055950A8	L070095A8	H080185B0.5
Size	Cover:105.5*15.5CM Heat patch:13*9.5CM	9.5*5.5CM	9.5*7CM	8*18.5CM
Weight	38±5G	30G	21G	11±2G
Material	Non-woven, PP Non-woven, Heating material	PP Non-woven, Heating material	PP Non-woven, Heating material, PE stick on clothes/PP film	Non-woven, Heating material
Using Position	Waist	Hand	Pelma	Eye

Item Number	H055095B8/L095130A12	L100150A8	H095285B8	H095285B8/H086260B8/H085250B8
Size	9.5*5.5CM/9.5*13CM	15*10CM	28.5*9.5CM	28.5*9.5/26.5*8.6/25*8.5CM
Weight	30±5G/38±5G	30G	55±5G	55±5G/51±5G/46±5G
Material	Spunlace PE/PP Non-woven, Heating material, PE Stick on skin/Spunlace stick on skin/ PE stick on clothes/PP film			PP Non-woven, Heating material, Ivory board film, PE stick on clothes
Using Position	Multi-use	Belly	Shoulder/Back	Foot



# 07 HOME CARE SERIES MORE CONVENIENT LIFE





## COOLING PATCH



Soothe fever either for Baby or Adult.  
Applicable to physically reduce the local temperature and prevent fever heat.  
Keep cool when driving or doing outdoor activities in summer.

### Advantages

Gentle to skin, strong adhesion, disposable.  
Soothe fever, cooling effect can last up to 8 hours.

		
Item Number	FCP 0305	FCP 0405
Size	12*5CM	13*9.5CM
Weight	10 G	18±2 G
Material	Non-woven	Non-woven







# ICE BAG

### Characteristic

Excellent cold therapy effects.  
Pleated design allows bag fabric to conform closely to body contours.  
Different sizes and colors available.  
Portable and reusable.

### Usage

Use with ice cubes, crushed ice or cold water.  
Inject the ice or cold water into the bag and close the bag tightly, then apply it on the body.



Item Number	ICE-BAG 5	ICE-BAG 6	ICE-BAG 9	ICE-BAG 11
Size	5 inch	6 inch	9 inch	11 inch
Weight	60G	68G	97G	120G
Material	polyester	polyester	polyester	polyester



# HOT WATER BAG

### Characteristic

Available in Plain and ribbed on one side.  
The hot water bag is made of natural rubber which is eco-friendly material, which is traditional product in winter season.  
Hot water bag can be used alone as a remedy for aches, pains and sports injuries.  
The bottle can also be used for cold purposes, allow easy filling of cold water and crushed ice.  
We can match different lovely or washable cotton/fleece/knitted cover, preventing the hot rubber from being in direct contact with skin.



Item Number	Size	Weight	Material
HWB001--2L	20*4*29CM	390G	Non-woven







# COOL GEL MAT

### Instruction

For pillow pad, just put it on your pillow. For bed mat, you can lie on it directly or put the bed sheet on the gel mat, then lie on it.

- Helps relieve night sweats and hot flashes
- Non-toxic and environmentally friendly
- Adopts special techniques, non-leaking and safe
- Enjoy coolness and softness at the same time
- Keep cool automatically, no electricity required, economical

Item Number	IBD 3040A	IBD 4343A	IBD 2740	IBD 4343B
Size	30*40CM	43*43CM	40*27CM	43*43CM
Weight	680G	810G	650G	650G
Material	Cotton/Polyester taffeta/Fabric	Fabric	Aluminum+Sponge	PVC



# COOLING GEL PILLOW

### Advantages

Lower your body temperature immediately.  
The cooling effect can lasts up to 3 hours.  
Non-toxic, Environmental-friendly, Energy-saving.  
Recent sleep studies have shown that cooling the brain is a simple, effective method to help relieve insomnia.  
Soothes sore, aching muscles, numbness, headaches and stress pain.  
Also great for hot flashes and nigh sweaters.  
Absorb heat from the body to help you sleep better in hot days.  
Ideal choice for home beds, toddler beds, hospital settings lounging around.

Item Number	CP 1040A	CP 1000	CUS 3035
Size	21.6*31.8CM	18.5*31CM	40*45CM
Weight	1040G	1000G	890G
Material	PVC+PE&NYLON	PE&NYLON+EVA	Polyester taffeta+Sponge

Item Number	CP 0812	CP 13201	CP 13202
Size	8*12CM	13*20CM	13*20CM
Weight	75G	270G	280G
Material	Polyester taffeta/Fabric	EVA	PVC compound spunlace

Reusable , keep in freezer for less than 2 hours before use.



# 08 MOM & BABY SERIES

Non-toxic, phthalate free, BPA free, latex free. Soft and flexible for comfort and feature a contoured shape that molds to your breast. With extra soft covers conform to breast and can be worn under bra.



## BREAST GEL PACK



### Advantages

Hot and cold therapy when Breast-feeding or during menstruation.

**Hot Therapy** to help relieve mastitis and help loosen clogged milk ducts to make feeding easier.

**Cold Therapy** to relieve pain, swelling, engorgement, plugged ducts, or mastitis.


### Breast Pump Aid

Use hot to reduce time spent pumping.


		
Item Number	BRA-GEL	BRA-GB
Size	17*17CM	17*17CM
Weight	165G	165G
Material	PVC+GEL	PVC+GEL Beads



Heat in water



Heat in microwave oven



Freeze in refrigerator





Advantages

- Perineal Instant Cold Pack especially used for gentle, effective postpartum care.
- Maximum comfort and protection with minimal nursing intervention.
- Cold therapy and absorbent pad in one.
- The all-in-one design eliminates the need for additional supplies, offering superior cost efficiency.
- Made of soft material, our perineal packs provide safe, controlled and effective thermal therapy and increased absorbency.

Alternate Activation

1. Position hands on inner water bag.
2. Firmly but gently twist inner bag to rupture.
3. Shake to mix contents.
4. Remove outer packaging and peel back adhesive liner.
5. Apply non-adhesive side to patient and secure pack to undergarment.



Item Number	Size	Weight	Material
PER 0613	13.5*34CM	220G	PE



Advantages

Strictly temperature control , gentle to baby skin.  
Disposable, Safe and no cross-infection.

Instruction

Bend metal disc back and forth with thumb and forefinger of both hands.  
Knead pack and secure warmer to heel by wrapping strip around foot.



Item Number	Size	Weight	Material
IHW 0305	12*8.5CM	57G	PE&Nylon












## 09 PROMOTION & GIFT SERIES

It's very convenient to use hand warmer outdoors, such as hiking. Bring you warmth in the cold weather. It can be made into various shape with different logo to meet your special needs. It's a good choice as gifts and promotional items with cheap price.

## CARTOON HOT/COLD TOY PACK

					
Item Number	HC-GEL 80	HC-GB 59	HC-GB 68	HC-GB 85	HC-GB 92
Size	10.7*10.7CM	10*10CM	11.6*9.7CM	11.9*10.5CM	7.6*11.6CM
Weight	80G	59G	68G	85G	92G
Material	PVC	PVC	PVC	PVC	PVC

					
Item Number	HC-GEL 88	HC-GEL 85	HC-GB 79	HC-GB 53	HC-GB 80
Size	9.8*10.8CM	9.7*10.6CM	10.3*10.3CM	7.6*10.5CM	9.5*10.8CM
Weight	88G	85G	79G	53G	80G
Material	PVC	PVC	PVC	PVC	PVC

					
Item Number	HC-GEL 122	HC-GEL 120	HC-GEL 135	HC-GEL 77	HC-GEL 138
Size	13*13CM	11*14.2CM	11*14.2CM	8*12CM	11.3*11.6CM
Weight	122G	120G	135G	77G	138G
Material	PVC	PVC	PVC	PVC	PVC